

## Product datasheet

### Anti-MLF2 antibody ab56378

#### 2 图像

#### 概述

产品名称	Anti-MLF2抗体
描述	小鼠单克隆抗体to MLF2
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein, corresponding to amino acids 1-249 of Human MLF2

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab56378** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 29 kDa.

#### 靶标

组织特异性	Ubiquitously expressed.
序列相似性	Belongs to the MLF family.
细胞定位	Cytoplasm. Nucleus.

## 图片



Anti-MLF2 antibody (ab56378) at 5 µg/ml +  
0.2 µg/lane of full length recombinant MLF2  
protein with a tag


### Secondary

Goat anti mouse IgG (H&L)-HRP conjugate

at 1/5000 dilution

**Predicted band size:** 29 kDa

**Observed band size:** 53 kDa



**All lanes :** Anti-MLF2 antibody (ab56378) at  
5 µg/ml

**Lane 1 :** 50 µg/lane of transfected 293T cell  
line lysate

**Lane 2 :** 50 µg/lane of non-transfected 293T  
cell line lysate

### Secondary

**All lanes :** Goat anti mouse IgG (H&L)-HRP  
conjugate at 1/2500 dilution

**Predicted band size:** 29 kDa

**Observed band size:** 30 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors