

Product datasheet

Anti-GRAP2 antibody ab56498

2 图像

概述

|       |   |
|-------|---|
| 产品名称  | Anti-GRAP2抗体  |
| 描述    | 小鼠单克隆抗体to GRAP2   |
| 宿主    | Mouse   |
| 经测试应用 | 适用于: WB, IHC-P  |
| 种属反应性 | 与反应: Human  |
| 免疫原   | Recombinant fragment: HHFHQERRGG SLDINDGHCG TGLGSEMNAALMHRRHDPVQLQAAGRVRWARALYDFEAL EDDELGFHSG EVVEVLDSSN PSWWTGRLHN , corresponding to amino acids 226-316 of Human GRAP2<br><a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a> |

性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None<br>PBS, pH 7.2   |
| 纯度   | Protein G purified  |
| 克隆   | 单克隆   |
| 同种型  | IgG2a   |
| 轻链类型 | kappa   |

应用

Our [Abpromise guarantee](#) covers the use of **ab56498** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| 应用 | Ab评论 | 说明  |
|----|------|---|
| WB |      | Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 38 kDa. |

| 应用    | Ab评论 | 说明                              |
|-------|------|---------------------------------|
| IHC-P |      | Use a concentration of 3 µg/ml. |

## 靶标


|              |  |
|--------------|--|
| <b>功能</b>    | Interacts with SLP-76 to regulate NF-AT activation. Binds to tyrosine-phosphorylated shc.  |
| <b>序列相似性</b> | Belongs to the GRB2/sem-5/DRK family.<br>Contains 1 SH2 domain.<br>Contains 2 SH3 domains. |

## 图片



GRAP2 antibody (ab56498) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human spleen.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-GRAP2 antibody (ab56498)



GRAP2 antibody (ab56498) at 1ug/lane + Jurkat cell lysate at 25ug/lane.

Western blot - Anti-GRAP2 antibody (ab56498)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors