

Product datasheet

Anti-LBX1 antibody ab56840

1 图像

概述

<b>产品名称</b>	Anti-LBX1抗体
<b>描述</b>	小鼠单克隆抗体to LBX1
<b>宿主</b>	Mouse
<b>经测试应用</b>	<b>适用于:</b> WB
<b>种属反应性</b>	<b>与反应:</b> Human
<b>免疫原</b>	Recombinant full length protein: 0 , corresponding to amino acids 133-219 of Human LBX1
	<a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>存储溶液</b>	Preservative: None PBS, pH 7.2
<b>纯度</b>	Protein G purified
<b>克隆</b>	单克隆
<b>同种型</b>	IgG2a
<b>轻链类型</b>	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab56840** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明


应用	Ab评论	说明
----	------	----

WB Use a concentration of 1 - 5 µg/ml.  
This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

## 靶标

细胞定位 Nuclear

## 图片



Western blot against tagged recombinant protein immunogen using ab56840 LBX1 antibody at 1ug/ml. Predicted band size of immunogen is 36 kDa

Western blot - Anti-LBX1 antibody (ab56840)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors