

Product datasheet

Anti-RAX antibody ab57183

3 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-RAX抗体 |
| 描述 | 小鼠单克隆抗体to RAX |
| 宿主 | Mouse |
| 经测试应用 | 适用于: WB, IHC-P |
| 种属反应性 | 与反应: Human |
| 免疫原 | Recombinant full length protein corresponding to Human RAX aa 1 0 4-2 0 6. Sequence: 0 Database link: Q9Y2V3 |

 [Run BLAST with](#)

 [Run BLAST with](#)

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None PBS, pH 7.2 |
| 纯度 | Protein G purified |
| 克隆 | 单克隆 |
| 同种型 | IgG2a |
| 轻链类型 | kappa |

应用

Our [Abpromise guarantee](#) covers the use of **ab57183** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|-------|------|--|
| WB | | Use a concentration of 1 - 5 µg/ml. |
| IHC-P | | Use at an assay dependent concentration. |
| IF | | Use at an assay dependent concentration. |


靶标

| | |
|--------------|---|
| 功能 | Plays a critical role in eye formation by regulating the initial specification of retinal cells and/or their subsequent proliferation. Binds to the photoreceptor conserved element-I (PCE-1/Ret 1) in the photoreceptor cell-specific arrestin promoter. |
| 组织特异性 | Expressed in the developing eye and weakly expressed in the adult retina. |
| 疾病相关 | Defects in RAX are the cause of microphthalmia isolated type 3 (MCOP3) [MIM:611038]. Microphthalmia is a clinically heterogeneous disorder of eye formation, ranging from small size of a single eye to complete bilateral absence of ocular tissues. Ocular abnormalities like opacities of the cornea and lens, scarring of the retina and choroid, cataract and other abnormalities like cataract may also be present. |
| 序列相似性 | Belongs to the paired homeobox family. Bicoid subfamily. Contains 1 homeobox DNA-binding domain. |
| 细胞定位 | Nucleus. |

图片




Western blot against tagged recombinant protein immunogen using ab57183 RAX antibody at 1ug/ml. Predicted band size of immunogen is 37 kDa



Immunohistochemical analysis of formalin-fixed paraffin-embedded human adrenal gland using ab57183 at concentration 3 ug/ml

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-RAX antibody (ab57183)



Immunofluorescence analysis of HeLa cell using ab57183 at concentration 10 ug/ml

Immunofluorescence - Anti-RAX antibody (ab57183)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors