

Product datasheet

Anti-Protocadherin beta 3 antibody ab57598

1 图像

概述

产品名称	Anti-Protocadherin beta 3抗体
描述	小鼠单克隆抗体to Protocadherin beta 3
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant fragment: ASEEIRKTFQ LNPITGDMQL VKYLNFEAIN SYEVDIEAKD GGGLSGKSTV IVQVVDVNDN PPELTLSSVN SPIPENSGET VLAVFSVSDL DSGDNGRV, corresponding to amino acids 284-382 of Human Protocadherin beta 3 Run BLAST with ExPASy Run BLAST with NCBI

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab57598** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

靶标

功能	Potential calcium-dependent cell-adhesion protein. May be involved in the establishment and maintenance of specific neuronal connections in the brain.
序列相似性	Contains 6 cadherin domains.
细胞定位	Cell membrane.

图片



Western blot against tagged recombinant protein immunogen using ab57598
Protocadherin beta 3 antibody at 1ug/ml.
Predicted band size of immunogen is 37 kDa

Western blot - Anti-Protocadherin beta 3 antibody (ab57598)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors