

Product datasheet

Anti-RASEF antibody ab58046

2 图像

概述

产品名称	Anti-RASEF 抗体
描述	小鼠单克隆抗体 to RASEF
宿主	Mouse
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant fragment, corresponding to amino acids 631-740 of Human RASEF

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG2b
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab58046** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 83 kDa.
IHC-P		Use a concentration of 3 µg/ml.

靶标


相关性	Binds predominantly GDP, and also GTP.
细胞定位	Cytoplasm, perinuclear region.

图片



RASEF antibody (ab58046) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human adrenal gland.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-RASEF antibody (ab58046)



RASEF antibody (ab58046) at 1ug/lane + K-562 cell lysate at 25ug/lane.

Western blot - Anti-RASEF antibody (ab58046)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors