

Product datasheet

Anti-CD90 / Thy1 antibody [30-H12] (FITC) ab62009

概述

产品名称	Anti-CD90 / Thy1 抗体[30-H12] (FITC)
描述	大鼠单克隆抗体[30-H12] to CD90 / Thy1 (FITC)
宿主	Rat
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: WB, IP, Flow Cyt, IHC-Fr
种属反应性	与反应: Mouse
免疫原	Full length mouse CD90/Thy1 protein.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 40% Glycerol, 1% BSA, PBS, pH 7.2
纯度	Protein G purified
纯化说明	ab62009 was purified from tissue culture supernatant.
克隆	单克隆
克隆编号	30-H12
同种型	IgG2b

应用

Our [Abpromise guarantee](#) covers the use of **ab62009** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
WB		Use at an assay dependent concentration. Predicted molecular weight: 18 kDa. This product is also suitable for use in identification of mouse T-cells subset population and phenotype analysis as a marker for progenitor cells during developmental stage.

