


Product datasheet

Anti-Apoptosis repressor with CARD antibody ab62149

3 图像

概述

产品名称	Anti-Apoptosis repressor with CARD抗体
描述	兔多克隆抗体 to Apoptosis repressor with CARD
宿主	Rabbit
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rhesus monkey 
免疫原	Synthetic peptide: SPSPTSRKGTSPKIPEGQ , corresponding to amino acids 191-208 of Human Apoptosis repressor with CARD Run BLAST with Run BLAST with
阳性对照	HeLa whole cell lysate. KB whole cell lysate. A549 whole cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab62149** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.
WB		1/1000.

靶标

功能	Isoform 1 may be involved in RNA splicing. Isoform 2 may inhibit apoptosis.
组织特异性	Highly expressed in heart and skeletal muscle. Detected at low levels in placenta, liver, kidney and pancreas.
序列相似性	Contains 1 CARD domain.
细胞定位	Cytoplasm and Nucleus > nucleolus.

图片



Western blot - Anti-Apoptosis repressor with CARD antibody (ab62149)

All lanes : Anti-Apoptosis repressor with CARD antibody (ab62149) at 1/1000 dilution


Lane 1 : HeLa whole cell lysate

Lane 2 : KB whole cell lysate

Lane 3 : A549 whole cell lysate


Observed band size: 25 kDa

Additional bands at: 28 kDa. We are unsure as to the identity of these extra bands.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Apoptosis repressor with CARD antibody (ab62149)

ab62149 at 181 ug/ml, staining Human Apoptosis repressor with CARD in heart tissue by IHC-P.



ab62149 at 5 μ g/ml, staining Human Apoptosis repressor with CARD in skeletal muscle by IHC-P.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Apoptosis repressor with CARD antibody (ab62149)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors