

Product datasheet

Anti-Phospholipase A2 activator protein (PLAP) antibody
ab63170

1 图像

概述

产品名称	Anti-Phospholipase A2活化剂蛋白(PLAP)抗体
描述	兔多克隆抗体to Phospholipase A2活化剂蛋白(PLAP)
宿主	Rabbit
经测试应用	适用于: ELISA, IHC-P
种属反应性	与反应: Human
免疫原	Synthetic peptide from the C-terminal region of human Phospholipase A2 activator protein (PLAP), conjugated to KLH.
阳性对照	Human lung carcinoma tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Protein A purified
纯化说明	ab63170 was purified through a protein A column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab63170** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/1000.
IHC-P		1/10 - 1/50.

靶标

功能	Involved in the maintenance of ubiquitin levels.
序列相似性	<p>Belongs to the WD repeat PLAP family.</p> <p>Contains 6 ARM repeats.</p> <p>Contains 1 PFU domain.</p> <p>Contains 1 PUL domain.</p> <p>Contains 7 WD repeats.</p>
结构域	<p>The PUL domain is composed of 6 armadillo-like repeats and mediates the interaction with VCP C-terminus.</p> <p>The PFU domain mediates interaction with ubiquitin.</p>

图片



ab63170, at 1/10 dilution, staining Phospholipase A2 activator protein (PLAP) in formalin-fixed, paraffin-embedded human lung carcinoma tissue by immunohistochemistry. ab63170 was peroxidase-conjugated to the secondary antibody, followed by DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Phospholipase A2 activator protein (PLAP) antibody (ab63170)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors