

Product datasheet

Anti-ZNF688 antibody ab66388

1 图像

概述

产品名称	Anti-ZNF688抗体
描述	兔多克隆抗体to ZNF688
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Synthetic peptide: KEEEPVPR AKGPRKAPVK ESPEVLVERN PDPAISVAPA RAQPPKNAAW, corresponding to a region within internal sequence amino acids 109-158 of Human ZNF688 <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	Jurkat cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab66388** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab 评论	说明

应用	Ab 评论	说明
WB		Use a concentration of 2.5 µg/ml. Detects a band of approximately 33 kDa (predicted molecular weight: 31 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
ELISA		Use at an assay dependent concentration. Titre using a peptide based assay: 1/312500.

## 靶标

功能	May be involved in transcriptional regulation.
序列相似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 2 C2H2-type zinc fingers. Contains 1 KRAB domain.
细胞定位	Nucleus.

## 图片



Anti-ZNF688 antibody (ab66388) at 2.5 µg/ml  
+ Jurkat cell lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 31 kDa

**Observed band size:** 33 kDa

**Additional bands at:** 15 kDa, 20 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors