

## Product datasheet

# Anti-MPPED1 antibody ab66988

### 1 图像

#### 概述

产品名称	Anti-MPPED1抗体
描述	小鼠多克隆抗体to MPPED1
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length protein, corresponding to amino acids 1-327 of Human MPPED1
阳性对照	Transfected 293T cell lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.5
纯度	Protein A purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab66988** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 37 kDa (predicted molecular weight: 37 kDa).

#### 靶标

功能	May have metallophosphoesterase activity (in vitro).
----	--


组织特异性

Expressed predominantly in adult brain.

序列相似性

Belongs to the UPF0046 family.

## 图片



**All lanes** : Anti-MPPED1 antibody (ab66988)  
at 1/500 dilution

**Lane 1** : MPPED transfected 293T cell lysate

**Lane 2** : non transfected 293T cell lysate

Lysates/proteins at 25 µg/ml per lane.

### Secondary

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP  
Conjugate secondary antibody at 1/2500  
dilution

**Predicted band size:** 37 kDa

**Observed band size:** 37 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors