

## Product datasheet

# Anti-MRPL40 antibody ab67184

### 1 图像

#### 概述

产品名称	Anti-MRPL40抗体
描述	小鼠多克隆抗体to MRPL40
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length MRPL40 protein (Human)
阳性对照	MRPL40 transfected 293T cell lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab67184** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

应用说明	WB: 1/500 - 1/1000. Detects a band of approximately 25 kDa (predicted molecular weight: 24 kDa). Detection Antibody Pair.
------	--

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

## 靶标


组织特异性

Ubiquitous.

细胞定位

Mitochondrion.

## 图片



**All lanes :** Anti-MRPL40 antibody (ab67184)

at 1/500 dilution

**Lane 1 :** MRPL40 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg/ml per lane.

### Secondary

**All lanes :** Goat Anti-Mouse IgG (H&L)-HRP

Conjugate at 1/2500 dilution

**Predicted band size:** 24 kDa

**Observed band size:** 25 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors