

Product datasheet

Anti-EFS antibody ab67502

★★★★☆ 1 Abreviews 1 图像

概述

产品名称	Anti-EFS抗体
描述	小鼠多克隆抗体to EFS
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Human
免疫原	Full length protein, corresponding to amino acids 1-468 of Human EFS
阳性对照	Transfected mammalian cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab67502** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB	★★★★☆	1/500 - 1/1000. Predicted molecular weight: 59 kDa.

靶标

相关性	EFS is a docking protein which plays a pivotal role for cell adhesion-related intracellular tyrosine kinase signaling, where it may serve as an activator of SRC and is known to interact with the SH3 domain of FYN and with CRK, SRC, and YES.
细胞定位	Cytoplasm

图片



Western blot - Anti-EFS antibody (ab67502)

All lanes : Anti-EFS antibody (ab67502) at 1/500 dilution

Lane 1 : EFS transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

Lane 1 : Goat Anti-Mouse IgG (H&L)-HRP Conjugated at 1/2500 dilution

Lane 2 : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

Predicted band size: 59 kDa

Observed band size: 70 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors