

Product datasheet

Anti-USP24 antibody ab67579

1 图像

概述

|       |                                    |
|-------|------------------------------------|
| 产品名称  | Anti-USP24抗体                       |
| 描述    | 小鼠多克隆抗体to USP24                    |
| 宿主    | Mouse                              |
| 经测试应用 | 适用于: WB                            |
| 种属反应性 | 与反应: Human                         |
| 免疫原   | Full length USP24 protein (Human)  |
| 阳性对照  | USP24 transfected 293T cell lysate |

性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None<br>Constituents: 1X PBS, pH 7.2  |
| 纯度   | Protein A purified  |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

应用

Our [Abpromise guarantee](#) covers the use of **ab67579** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明  |
|----|------|---|
| WB |      | 1/500 - 1/1000. Detects a band of approximately 37 kDa (predicted molecular weight: 211 kDa). |

靶标

|              |   |
|--------------|---|
| <b>功能</b>    | Involves in the ubiquitin-dependent proteolytic pathway in conjunction with the 26S proteasome. |
| <b>序列相似性</b> | Belongs to the peptidase C19 family.<br>Contains 1 UBA domain.                                  |
| <b>翻译后修饰</b> | Phosphorylated upon DNA damage, probably by ATM or ATR.   |

## 图片



**All lanes :** Anti-USP24 antibody (ab67579) at 1/500 dilution

**Lane 1 :** USP24 transfected 293T cell lysate

**Lane 2 :** Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Predicted band size:** 211 kDa

**Observed band size:** 37 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors