

Product datasheet

Anti-DCAKD antibody ab67667

1 图像

概述

产品名称	Anti-DCAKD抗体
描述	小鼠多克隆抗体to DCAKD
宿主	Mouse
特异性	This antibody is reactive against mammalian transfected lysate.
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length DCAKD protein (Human)
阳性对照	DCAKD transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab67667** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

**应用说明** WB: 1/500 - 1/1000. Detects a band of approximately 24 kDa. Detection Antibody Pair.

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

## 靶标

### 序列相似性

Belongs to the CoaE family.  
Contains 1 DPCK (dephospho-CoA kinase) domain.

## 图片



**All lanes :** Anti-DCAKD antibody (ab67667) at 1/500 dilution

**Lane 1 :** DCAKD transfected 293T cell lysate

**Lane 2 :** Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Observed band size:** 24 kDa

[why is the actual band size different from the predicted?](#)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors