

Product datasheet

Anti-LYPD5 antibody ab68079

1 图像

概述

|       |  |
|-------|--|
| 产品名称  | Anti-LYPD5抗体   |
| 描述    | 小鼠多克隆抗体to LYPD5  |
| 宿主    | Mouse  |
| 经测试应用 | 适用于: WB  |
| 种属反应性 | 与反应: Human   |
| 免疫原   | Full length protein, corresponding to amino acids 1-208 of Human LYPD5 |
| 阳性对照  | Transfected 293T cell lysate   |

性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None<br>Constituents: PBS, pH 7.2   |
| 纯度   | Protein A purified  |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

应用

Our [Abpromise guarantee](#) covers the use of **ab68079** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明   |
|----|------|--|
| WB |      | 1/500 - 1/1000. Detects a band of approximately 23 kDa (predicted molecular weight: 27 kDa). |


靶标

**序列相似性**

Contains 1 UPAR/Ly6 domain.

**细胞定位**

Cell membrane.

**图片**

Western blot - Anti-LYPD5 antibody (ab68079)

**All lanes :** Anti-LYPD5 antibody (ab68079) at 1/500 dilution

**Lane 1 :** LYPD5 transfected 293T cell lysate

**Lane 2 :** non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

**Secondary**

**All lanes :**

Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/25 dilution

**Predicted band size:** 27 kDa

**Observed band size:** 23 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors