

Product datasheet

Anti-PATZ1 antibody ab68646

1 图像

概述

产品名称	Anti-PATZ1抗体
描述	小鼠多克隆抗体to PATZ1
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	PATZ1 full length (NP_114440.1) human protein
阳性对照	PATZ1 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.4 Constituent: PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab68646** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 67 kDa (predicted molecular weight: 74 kDa).

靶标

<b>功能</b>	Transcriptional repressor.
<b>组织特异性</b>	Ubiquitous.
<b>疾病相关</b>	A chromosomal aberration involving PATZ1 is associated with small round cell sarcoma. Translocation t(1;22)(p36.1;q12) with EWSR1.
<b>序列相似性</b>	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 1 A.T hook DNA-binding domain. Contains 1 BTB (POZ) domain. Contains 7 C2H2-type zinc fingers.
<b>细胞定位</b>	Nucleus.

## 图片



Western blot - Anti-PATZ1 antibody (ab68646)

**All lanes** : Anti-PATZ1 antibody (ab68646) at 1/500 dilution

**Lane 1** : PATZ1 transfected 293T cell lysate

**Lane 2** : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP Conjugate at 1/2500 dilution

**Predicted band size:** 74 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors