

Product datasheet

Anti-THAP4 antibody ab68744

1 图像

概述

产品名称	Anti-THAP4抗体
描述	小鼠多克隆抗体to THAP4
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length protein, corresponding to amino acids 1-577 of Human THAP4
阳性对照	THAP4 transfected 293T cell line

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab68744** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 63 kDa.

靶标


序列相似性

Contains 1 THAP-type zinc finger.

翻译后修饰

Phosphorylated upon DNA damage, probably by ATM or ATR.

图片



Western blot - Anti-THAP4 antibody (ab68744)

All lanes : Anti-THAP4 antibody (ab68744) at 1/500 dilution

Lane 1 : transfected 293T cell line

Lane 2 : non transfected 293T cell line

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP

Conjugate secondary antibody at 1/2500 dilution

Predicted band size: 63 kDa

Observed band size: 75 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors