

Product datasheet

Anti-DFFB/CAD antibody [4B8] ab71083

1 图像

概述

产品名称	Anti-DFFB/CAD抗体[4B8]
描述	小鼠单克隆抗体[4B8] to DFFB/CAD
宿主	Mouse
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 不与反应: Mouse, Rat
免疫原	Recombinant full length human DFFB purified from E.coli with His tag.
阳性对照	A431, Jurkat, HeLa and HepG2 cell lysates.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	Preservative: 0.03% Sodium azide Constituents: HEPES, 0.01% BSA, 50% Glycerol, 0.87% Sodium chloride
纯度	Ammonium Sulphate Precipitation
克隆	单克隆
克隆编号	4B8
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab71083** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000. Detects a band of approximately 39 kDa (predicted molecular weight: 39 kDa).
ELISA		Use at an assay dependent dilution.

靶标

功能	Nuclease that induces DNA fragmentation and chromatin condensation during apoptosis. Degrades naked DNA and induces apoptotic morphology.
序列相似性	Contains 1 CIDE-N domain.
细胞定位	Cytoplasm. Nucleus.

图片



All lanes : Anti-DFFB/CAD antibody [4B8] (ab71083) at 1/1000 dilution

Lane 1 : A431 cell lysate

Lane 2 : Jurkat T cell lysate

Lane 3 : HeLa cell lysate

Lane 4 : HepG2 cell lysate

Predicted band size: 39 kDa

Observed band size: 39 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors