

Product datasheet

Anti-ART3 antibody ab71638

1 图像

概述

产品名称	Anti-ART3抗体
描述	兔多克隆抗体to ART3
宿主	Rabbit
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide from the N-terminal region of human ART3, conjugated to KLH
阳性对照	Mouse brain tissue lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Protein G purified
纯化说明	The antibody was purified through a protein G column, eluted with high and low pH buffers and neutralized immediately, followed by dialysis against PBS.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab71638** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/1000.
WB		1/100 - 1/500. Predicted molecular weight: 44 kDa.

## 靶标

组织特异性	Testis specific.
序列相似性	Belongs to the Arg-specific ADP-ribosyltransferase family.
细胞定位	Cell membrane.

## 图片



**All lanes :** Anti-ART3 antibody (ab71638) at 1/100 dilution

**Lane 1 :** Mouse brain tissue lysate

**Lane 2 :** Mouse brain tissue lysate

Lysates/proteins at 12.5 µg per lane.

**Predicted band size:** 44 kDa

**Observed band size:** 44 kDa

**Additional bands at:** 100 kDa, 150 kDa, 28 kDa, 32 kDa, 50 kDa, 55 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors