

Product datasheet

Goat Anti-Mouse IgG1 (SureLight® Allophycocyanin) ab72553

概述

产品名称	山羊抗小鼠IgG1 (SureLight® Allophycocyanin荧光)
宿主	Goat
靶标种属	Mouse
经测试应用	适用于: Flow Cyt
免疫原	Mouse IgG1
偶联物	SureLight® Allophycocyanin. Ex: 652nm, Em: 657nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 2mM Sodium Azide Constituents: 0.1% BSA, 50mM Sucrose, 150mM Sodium chloride, 100mM Sodium phosphate, pH 7.4
纯化说明	Size exclusion chromatography as well as some proprietary methods for purifying the final materials were used.
共轭说明	Excitation max. λ 650 nm. Emission max. λ 660 nm.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab72553** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors