

Product datasheet

Anti-PABPC5 antibody ab72861

1 图像

概述

产品名称	Anti-PABPC5抗体
描述	小鼠多克隆抗体to PABPC5
宿主	Mouse
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Recombinant full length Humanprotein (NP_543022)
阳性对照	Reactive against mammalian transfected lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab72861** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

ELISA

应用说明	ELISA: Use at an assay dependent dilution. WB: 1/500 - 1/1000. Predicted molecular weight: 43 kDa. This antibody has only been tested in
------	---

WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

Not yet tested in other applications.


Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性 PABPC5 (polyadenylate-binding protein 5) binds the polyadenylate tail of mRNA and is thought to be involved in cytoplasmic regulatory processes of mRNA metabolism. In vivo, it may also bind to cytoplasmic RNA sequences other than polyadenylate.

细胞定位 Cytoplasmic

图片



All lanes : Anti-PABPC5 antibody (ab72861)
at 1/500 dilution

Lane 1 : PABPC5-transfected 293T cell
lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L) HRP-
conjugated at 1/5000 dilution

Predicted band size: 43 kDa

Observed band size: 43 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors