

Product datasheet

Anti-ZNHIT3 antibody ab74195

1 图像

概述

产品名称	Anti-ZNHIT3抗体
描述	兔多克隆抗体to ZNHIT3
宿主	Rabbit
经测试应用	适用于: IP 不适用于: WB
种属反应性	与反应: Human 预测可用于: Pig, Chimpanzee
免疫原	Synthetic peptide corresponding to a region between residues 25 and 75 of human ZNHIT3 (NP_004764.1)
阳性对照	HeLa whole cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 8mM PBS, 60mM Citrate, 150mM Tris, pH 7-8
纯度	Immunogen affinity purified
纯化说明	ab74195 was affinity purified using an epitope specific to ZNHIT3 immobilized on solid support.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab74195** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IP		Use at 2-5 µg/mg of lysate.

应用说明


Is unsuitable for WB.

靶标

相关性

ZNHIT3 functions as a transcriptional coactivator and contains 1 HIT type zinc finger.

图片



Detection of ZNHIT3 by Western Blot of Immunoprecipitate.
anti-ZNHIT3 (ab74196), at 1 µg/ml, staining ZNHIT3 in HeLa whole cell lysate immunoprecipitated using ab74195 at 3 µg/mg lysate (1 mg/IP; 20% of IP loaded/lane).
Detection: Chemiluminescence with exposure time of 30 seconds.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors