

Product datasheet

Anti-ALS2CR4 antibody ab74236

★★★★☆ 1 Abreviews 1 图像

概述

产品名称	Anti-ALS2CR4抗体
描述	兔多克隆抗体to ALS2CR4
宿主	Rabbit
特异性	This antibody detects endogenous levels of total ALS2CR4 protein.
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Mouse, Human
免疫原	A synthesized peptide derived from the internal region of human ALS2CR4.
阳性对照	Extracts from Jurkat cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab74236** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB	★★★★☆	1/500 - 1/1000. Detects a band of approximately 48 kDa (predicted molecular weight: 48 kDa).
ELISA		1/40000.

靶标


相关性

There are two named isoforms of ALS2CR4.

细胞定位

Membrane; Multi-pass membrane protein

图片



All lanes : Anti-ALS2CR4 antibody (ab74236) at 1/500 dilution

Lane 1 : Extracts from Jurkat cells

Lane 2 : Extracts from Jurkat cells plus 5µg immunizing peptide

Lysates/proteins at 5 µg per lane.

Predicted band size: 48 kDa

Observed band size: 48 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors