


Product datasheet

Anti-FUBP3 antibody ab74570

1 图像

概述

产品名称	Anti-FUBP3抗体
描述	兔多克隆抗体to FUBP3
宿主	Rabbit
特异性	ab74570 detects endogenous levels of total FUBP3 protein.
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthetic peptide derived from internal sequence of human FUBP3.
阳性对照	Extracts from COLO cells and HepG2 cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg ²⁺ and Ca ²⁺), 150mM Sodium chloride, pH 7.4
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab74570** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 61 kDa.
ELISA		1/20000.

靶标

功能	May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.
组织特异性	Detected in a number of cell lines.
序列相似性	Contains 4 KH domains.
翻译后修饰	Phosphorylated upon DNA damage, probably by ATM or ATR.
细胞定位	Nucleus.

图片



Western blot - Anti-FUBP3 antibody (ab74570)

All lanes : Anti-FUBP3 antibody (ab74570) at 1/500 dilution

Lane 1 : Extracts from COLO cells

Lane 2 : Extracts from HepG2 cells

Lane 3 : Extracts from COLO cells with immunising peptide at 10 μ g

Lysates/proteins at 30 μ g per lane.

Predicted band size: 61 kDa

Observed band size: 61 kDa

Additional bands at: 28 kDa, 50 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors