


Product datasheet

Anti-KLKBL4 antibody ab76035

1 图像

概述

产品名称	Anti-KLKBL4抗体
描述	兔多克隆抗体to KLKBL4
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Horse 
免疫原	Synthetic peptide (Human)
阳性对照	Fetal muscle lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab76035** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 44 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

序列相似性	Belongs to the peptidase S1 family. Plasma kallikrein subfamily. Contains 1 peptidase S1 domain.
结构域	The serine protease domain lacks the essential His, Asp, and Ser residues of the catalytic triad at positions 83, 129 and 221 and is predicted to be inactive.
细胞定位	Secreted.

## 图片



Anti-KLKBL4 antibody (ab76035) at 1  $\mu$ g/ml  
(in 5% skim milk / PBS buffer) + fetal muscle  
lysate at 10  $\mu$ g

### Secondary

HRP-conjugated anti-Rabbit IgG at 1/50000  
dilution

**Predicted band size:** 44 kDa

**Observed band size:** 40 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors