

Product datasheet

Anti-ZPLD1 antibody ab76704

1 图像

概述

产品名称	Anti-ZPLD1抗体
描述	小鼠单克隆抗体to ZPLD1
宿主	Mouse
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Recombinant fragment 预测可用于: Mouse, Rat, Horse, Cow, Dog, Human, Pig, Chimpanzee, Cynomolgus monkey, Macaque monkey 
免疫原	GST-tagged recombinant fragment: WNVLMDYCYT TPSGNPNDI RYDLFLSCDK DPQTTVIENG RSQRGRFSFE VFRFVKHKNQ KMSTVFLHCV TKLCRADD CP FLMPICSHRE RRDAGRRTTW, corresponding to amino acids 244-344 of Human ZPLD1. <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab76704** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

ELISA

#### 应用说明

ELISA: Use at an assay dependent dilution.

WB: 1/500 - 1/1000. Predicted molecular weight: 46 kDa. Note: the antibody was only tested against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

#### 靶标


##### 序列相似性

Contains 1 ZP domain.

##### 细胞定位

Membrane.

#### 图片



Anti-ZPLD1 antibody (ab76704) at 1/500 dilution + immunogen (GST-tagged recombinant fragment of human ZPLD1) at 0.2 µg

##### Secondary

Goat Anti-Mouse IgG (H&L)-HRP at 1/5000 dilution

**Predicted band size:** 46 kDa

MW of GST tag: 26 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors