

Product datasheet

Anti-POU5F2 antibody ab76705

1 图像

概述

产品名称	Anti-POU5F2抗体
描述	小鼠单克隆抗体to POU5F2
宿主	Mouse
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Recombinant fragment 预测可用于: Human, Chimpanzee
免疫原	GST-tagged recombinant fragment: AGHRPSNHFC PLPGSGGGGP RGPMLRVDT LTWLSTQAAP GRVMVWPAVR PGICPGPDVW RIPLGPLPHE FRGWIAPCRP RLGASEAGDW LRRPSEGAL, corresponding to amino acids 2-100 of Human FLJ25680 Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	购买相配的WB阳性对照 Recombinant Human POU5F2 protein >

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab76705** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		
ELISA		


应用说明

ELISA: Use at an assay dependent dilution.
 WB: 1/500 - 1/1000. Predicted molecular weight: 36 kDa. Note: the antibody was only tested against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

Not yet tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

靶标

图片



Western blot - Anti-POU5F2 antibody (ab76705)

Anti-POU5F2 antibody (ab76705) at 1/500 dilution + immunogen (GST-tagged recombinant fragment of Human FLJ25680) at 0.2 µg

Secondary

Goat Anti-Mouse IgG (H&L)-HRP at 1/5000 dilution

Predicted band size: 36 kDa

MW of GST tag: 26 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors