

Product datasheet

Anti-ADH1C antibody ab76745

2 图像

概述

产品名称	Anti-ADH1C抗体
描述	小鼠多克隆抗体to ADH1C
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length Human ADH1C protein (NP_000660.1, 1 a.a. ~ 375 a.a).
阳性对照	ADH1C transfected 293T cell lysate; human liver tissue lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab76745** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Detects a band of approximately 40 kDa (predicted molecular weight: 40 kDa).

靶标

序列相似性 Belongs to the zinc-containing alcohol dehydrogenase family.

图片



Anti-ADH1C antibody (ab76745) at 1 μ g/ml + Human liver tissue lysate at 25 μ g


Secondary

Goat Anti-Mouse IgG (H&L)-HRP at 1/2500 dilution

Predicted band size: 40 kDa

Observed band size: 37 kDa

[why is the actual band size different from the predicted?](#)



All lanes : Anti-ADH1C antibody (ab76745) at 1 μ g/ml

Lane 1 : ADH1C transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 μ g per lane.

Secondary

All lanes : Goat Anti-Mouse IgG (H&L)-HRP at 1/2500 dilution

Predicted band size: 40 kDa

Observed band size: 40 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors