

Product datasheet

Anti-RNF180 antibody ab76803

1 图像

概述

产品名称	Anti-RNF180抗体
描述	小鼠单克隆抗体to RNF180
宿主	Mouse
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Recombinant fragment 预测可用于: Human, Orangutan
免疫原	Recombinant fragment: MKRSKELITK NHSQEETSIL RCWKCRKCIA SSGCFMEYLE NQVIKDKDDS VDAQNICHVW HMNVEALPEW ISCLIQKAQW TVGKLNCPFC GARLGGFNFV , corresponding to amino acids 1-101 of Human RNF180(NP_848627), with a 26kDa enzyme tag. Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Recombinant fragment of human RNF180.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
同种型	IgG2b
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab76803** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		
ELISA		

应用说明

ELISA: Use at an assay dependent dilution.
 WB: Use at a concentration of 1 - 5 µg/ml. Detects a band of approximately 37 kDa (predicted molecular weight: 68 kDa).


This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

Not yet tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性 RNF180 (RING finger protein 180) contains 1 RING-type zinc finger.
细胞定位 Membrane; Single-pass membrane protein

图片



Western blot - Anti-RNF180 antibody (ab76803)

Anti-RNF180 antibody (ab76803) at 5 µg/ml +
 Recombinant fragment of human RNF180 at
 0.2 µg

Secondary

Goat Anti-Mouse IgG (H&L)-HRP at 1/5000
 dilution

Predicted band size: 68 kDa

Observed band size: 37 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors