

Product datasheet

Anti-LRRC24 antibody ab83257

1 图像

概述

产品名称	Anti-LRRC24抗体
描述	兔多克隆抗体to LRRC24
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Horse, Guinea pig, Cow, Cat, Dog
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 72-121 (PLAALRRLYL HNNSLRALEA GAFRAQPRLL ELALTSNRLR GLRSGAFVGL) of Human LRRC24 (NP_001019849) Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	HeLa cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab83257** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 55 kDa (predicted molecular weight: 55 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
ELISA		Use at an assay dependent concentration. Titre using peptide based assay: 1/312500.

靶标

相关性	LRRC24 contains 1 Ig-like C2-type (immunoglobulin-like) domain and 7 LRR (leucine-rich) repeats. It is a single-pass membrane protein. The function of the LRRC24 protein remains unknown.
细胞定位	Membrane; Single-pass membrane protein

图片



Anti-LRRC24 antibody (ab83257) at 1 µg/ml
(5% skim milk / PBS buffer) + HeLa cell lysate
at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 55 kDa

Observed band size: 55 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors