

Product datasheet

Anti-ZFP37 antibody ab84579

1 图像

概述

产品名称	Anti-ZFP37抗体
描述	兔多克隆抗体to ZFP37
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Cow, Cat, Dog, Pig, Chimpanzee
免疫原	Synthetic peptide corresponding to a region within internal sequence amino acids 252-301 (LCCHSSSHIK QDKIQTGEKH EKSPSLSSST KHEKPQACVK PYECNQCGKV) of Human ZFP37 (NP_003399). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	721_B cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab84579** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----


应用	Ab评论	说明
----	------	----

WB Use a concentration of 1 µg/ml. Predicted molecular weight: 71 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

相关性 ZFP37 may be involved in transcriptional regulation.
细胞定位 Nuclear

图片



Anti-ZFP37 antibody (ab84579) at 1 µg/ml (in 5% skim milk / PBS buffer) + 721_B cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 71 kDa

Observed band size: 68 kDa

Additional bands at: <22 kDa. We are unsure as to the identity of these extra bands.

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors