


Product datasheet

Anti-ANKRD13B antibody ab85577

1 图像

概述

产品名称	Anti-ANKRD13B抗体
描述	兔多克隆抗体to ANKRD13B
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Guinea pig, Cow, Cat, Dog 
免疫原	Synthetic peptide corresponding to a region within the middle amino acids 540-589 (HPMSYEGRRQDRSAPPTPQRQPAPPASVPSRPSSGPGSGGHVFRSYDEQ) of Human ANKRD13B (NP_689558) Run BLAST with Run BLAST with
阳性对照	Human fetal heart lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab85577** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 70 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

相关性

The exact function of ANKRD13B remains unknown.

图片



Anti-ANKRD13B antibody (ab85577) at 1 µg/ml (in 5% skim milk / PBS buffer) + Fetal heart lysate (human) at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 70 kDa

Observed band size: 70 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors