

Product datasheet

Anti-IgG antibody [9E2] ab85791

1 References

概述

产品名称	Anti-IgG抗体[9E2]
描述	小鼠单克隆抗体[9E2] to IgG
宿主	Mouse
特异性	Detects specific 'late' antibodies.
经测试应用	适用于: Indirect ELISA
种属反应性	与反应: Sheep, Goat, Cow, Deer
免疫原	Full length purified native protein from Sheep

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.4
纯度	Protein A purified
克隆	单克隆
克隆编号	9E2
骨髓瘤	Sp2/0
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab85791** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Indirect ELISA		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors