

Product datasheet

Anti-AP3M2 antibody ab87092

1 References 1 图像

概述

产品名称	Anti-AP3M2抗体
描述	兔多克隆抗体to AP3M2
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Pig, Zebrafish



免疫原	Synthetic peptide from within residues: VNTITGSTNV GDQLPTGQLS VVPWRRRTGVK YTNNEAYFDV IEEIDAIIDK, corresponding to internal sequence amino acids 144-193 of Human AP3M2 (NP_006794).
-----	---

 [Run BLAST with ExPASy](#)  [Run BLAST with NCBI](#)

阳性对照	HepG2 cell lysate
------	-------------------

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab87092** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 47 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

相关性 Part of the AP-3 complex, AP3M3 is an adaptor-related complex not clathrin-associated. The complex is associated with the Golgi region and to more peripheral structures. It facilitates the budding of vesicles from the Golgi membrane and may be directly involved in trafficking to lysosomes.

细胞定位 Golgi apparatus.

图片



Anti-AP3M2 antibody (ab87092) at 1 µg/ml +
HepG2 cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 47 kDa

Observed band size: 47 kDa

Additional bands at: 65 kDa. We are unsure
as to the identity of these extra bands.

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors