

Product datasheet

Anti-CLYBL antibody ab88571

2 图像

概述

产品名称	Anti-CLYBL抗体
描述	小鼠多克隆抗体to CLYBL
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length Human protein (AAH34360.1)
阳性对照	Liver lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab88571** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 37 kDa.

靶标

相关性 Belongs to the hpch/hpal aldolase family, citrate lyase beta subunit-like subfamily. Binds 1 magnesium ion per subunit. There are two isoforms.


细胞定位 Mitochondrial

图片



Anti-CLYBL antibody (ab88571) at 1 $\mu\text{g/ml}$ +
Human liver lysate at 50 μg

Predicted band size: 37 kDa



All lanes : Anti-CLYBL antibody (ab88571) at
1 $\mu\text{g/ml}$

Lane 1 : transfected 293T lysate

Lane 2 : non-transfected lysate

Lysates/proteins at 25 μg per lane.

Predicted band size: 37 kDa

Observed band size: 37 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors