

Product datasheet

Anti-SFXN2 antibody ab88916

2 图像

概述

产品名称	Anti-SFXN2抗体
描述	小鼠多克隆抗体to SFXN2
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length human SFXN2 protein (AAH22091).
阳性对照	Human kidney tissue lysate and SFXN2 transfected 293T cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab88916** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 36 kDa.

靶标


相关性

SFXN2 may be a potential iron transporter.

细胞定位

Mitochondrion membrane; Multi pass membrane


图片



Western blot - Anti-SFXN2 antibody (ab88916)

Anti-SFXN2 antibody (ab88916) at 1 µg/ml + human kidney tissue lysate at 50 µg

Predicted band size: 36 kDa



Western blot - Anti-SFXN2 antibody (ab88916)

All lanes : Anti-SFXN2 antibody (ab88916) at 1 µg/ml

Lane 1 : SFXN2 transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Predicted band size: 36 kDa

Observed band size: 41 kDa

[why is the actual band size different from the predicted?](#)

Additional bands at: 26 kDa, 28 kDa, 35 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors