

Product datasheet

Anti-C9orf68 antibody ab89115

2 图像

概述

产品名称	Anti-C9orf68抗体
描述	小鼠多克隆抗体to C9orf68
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length Human C9orf68 protein (AAH34293).
阳性对照	LNCaP cell lysate and C9orf68 transfected 293T cell line lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab89115** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 50 kDa (predicted molecular weight: 45 kDa).

靶标

## 相关性

The function of this protein is still undiscovered.


## 图片



Anti-C9orf68 antibody (ab89115) at 1  $\mu$ g/ml +  
LNCaP Cell lysate at 50  $\mu$ g

**Predicted band size:** 45 kDa

**Observed band size:** 50 kDa



**All lanes :** Anti-C9orf68 antibody (ab89115)  
at 1  $\mu$ g/ml

**Lane 1 :** C9orf68 transfected 293T cell line  
lysate

**Lane 2 :** Non-transfected lysate

Lysates/proteins at 25  $\mu$ g per lane.

**Predicted band size:** 45 kDa

**Observed band size:** 50 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors