

Product datasheet

Anti-SVH antibody ab89319

2 图像

概述

产品名称	Anti-SVH抗体
描述	小鼠多克隆抗体to SVH
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length Human SVH (AAH03586, amino acids 1-308)
阳性对照	Human liver tissue lysate and SVH transfected 293T cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab89319** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

应用说明 WB: 1/500 - 1/1000. Predicted molecular weight: Isoform 1, 38 kDa; Isoform 2, 34 kDa.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能	May play a role in cell survival and cell growth. May suppress the transcriptional activity of p53/TP53.
组织特异性	Expressed in all tissues tested with higher expression in placenta, liver, kidney, heart and brain.
序列相似性	Contains 1 ARM repeat.
细胞定位	Endoplasmic reticulum membrane.

图片




Anti-SVH antibody (ab89319) at 1/500 dilution + Human liver tissue lysate at 50 μ g

Predicted band size: 34 kDa

Observed band size: 30 kDa

[why is the actual band size different from the predicted?](#)



All lanes : Anti-SVH antibody (ab89319) at 1/500 dilution

Lane 1 : SVH transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 μ g per lane.

Predicted band size: 34 kDa

Observed band size: 28 kDa [why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors