

Product datasheet

Anti-KRR1 antibody ab89601

3 图像

概述

产品名称	Anti-KRR1抗体
描述	小鼠多克隆抗体to KRR1
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Full length Human KRR1 (AAH26107.1)
阳性对照	Human pancreas tissue lysate. HeLa cell lysate. KRR1 transfected 293T cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 1X PBS, pH 7.2
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab89601** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 44 kDa.

靶标

功能	Required for 40S ribosome biogenesis. Involved in nucleolar processing of pre-18S ribosomal RNA and ribosome assembly.
序列相似性	Belongs to the KRR1 family. Contains 1 KH domain.
细胞定位	Nucleus > nucleolus.

图片




Western blot - Anti-KRR1 antibody (ab89601)

Anti-KRR1 antibody (ab89601) at 1 µg/ml +
Human pancreas tissue lysate at 50 µg

Predicted band size: 44 kDa

Observed band size: 44 kDa




Western blot - Anti-KRR1 antibody (ab89601)

Anti-KRR1 antibody (ab89601) at 1 µg/ml +
HeLa cell lysate at 50 µg

Predicted band size: 44 kDa

Observed band size: 50 kDa



Western blot - Anti-KRR1 antibody (ab89601)

All lanes : Anti-KRR1 antibody (ab89601) at
1 µg/ml

Lane 1 : KRR1 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Predicted band size: 44 kDa

Observed band size: 48 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors