

Product datasheet

Anti-IL17A Receptor antibody [MM0371-8D15] ab92608

1 References

概述

产品名称	Anti-IL17A Receptor抗体[MM0371-8D15]
描述	小鼠单克隆抗体[MM0371-8D15] to IL17A Receptor
宿主	Mouse
经测试应用	适用于: WB, Neutralising
种属反应性	与反应: Human
免疫原	Recombinant fragment, corresponding to the extracellular domain of Human IL17A Receptor.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS
纯度	Protein G purified
纯化说明	Purified from tissue culture supernatant and 0.2 µm filtered.
克隆	单克隆
克隆编号	MM0371-8D15
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab92608** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

WB

Neutralising

应用说明 Neut: Use at an assay dependent dilution.

WB: 1/500 - 1/1000. Predicted molecular weight: 96 kDa.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

功能	Receptor for IL17A and IL17F. Binds its IL17A ligand with low affinity, suggesting that additional components are involved in IL17A-induced signaling.
组织特异性	Widely expressed.
疾病相关	Defects in IL17RA are the cause of familial candidiasis type 5 (CANDF5) [MIM:613953]. CANDF5 is a rare disorder with altered immune responses and impaired clearance of fungal infections, selective against Candida. It is characterized by persistent and/or recurrent infections of the skin, nails and mucous membranes caused by organisms of the genus Candida, mainly Candida albicans.
序列相似性	Contains 1 SEFIR domain.
翻译后修饰	Glycosylated.
细胞定位	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors