

Product datasheet

Anti-Znt6 antibody ab93249

概述

产品名称	Anti-Znt6抗体
描述	兔多克隆抗体to Znt6
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Recombinant fragment 预测可用于: Mouse, Rat
免疫原	Synthetic peptide corresponding to the very C terminal residues of Mouse Znt6 (Zinc transporter 6) protein.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.01% Sodium Azide Constituents: 50% Glycerol, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab93249** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
----	-------	----

WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 51 kDa. This antibody has been tested in Western blot against the recombinant peptide used as an immunogen. We have no data on detection of endogenous protein.
----	--	---

应用	Ab 评论	说明
----	----------	----

ELISA 1/2000 - 1/5000.

靶标

功能	Zinc-efflux transporter which allocates the cytoplasmic zinc to the trans-Golgi network (TGN) as well as the vesicular compartment.
组织特异性	Expressed in brain; especially in cerebellum, hippocampus, parahippocampal gyrus, superior and middle temporal gyrus. Also expressed in B cells, colon, eye, and lung. Lower expression was present in bone, brain, cervix, ear, heart, kidney, muscle, nerve, pancreas, prostate, skin, stomach, and testis.
序列相似性	Belongs to the cation diffusion facilitator (CDF) transporter (TC 2.A.4) family. SLC30A subfamily.
细胞定位	Golgi apparatus > trans-Golgi network membrane. Found in vesicles.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors