

## Product datasheet

# Antigen Retrieval Buffer (100X Citrate Buffer pH 6.0) ab94674

### 1 References

#### 概述

<b>产品名称</b>	Antigen Retrieval缓冲液(100X Citrate缓冲液pH 6.0)
<b>描述</b>	Antigen Retrieval缓冲液(100X Citrate缓冲液pH 6.0)
<b>经测试应用</b>	<b>适用于:</b> IHC-P
<b>常规说明</b>	<p><b>PROTOCOL:</b></p> <p>Place slides in a Coplin jar containing 1x Target Retrieval Buffer; cover with a vented plastic wrap and place the jar in microwave and set high power to boil and set low power to keep it boiling for 10 min. Let the sections cool in the microwave for at least 20min.</p> <p>Wash sections with hot tap water for 1 minute.</p> <p>Wash sections in buffer for 2x3 minutes.</p> <p>Block endogenous peroxidase as usual.</p> <p>Wash sections in buffer for 2x3 minutes.</p> <p>Block non-specific sites with a Blocking reagent as usual.</p> <p>Place optimally diluted primary antibody on the sections (incubation time and temperature for a given set of experimental conditions should be determined by the investigator).</p> <p>Wash sections in buffer for 2x3 minutes</p> <p>Rest of the procedure is same as routinely performed in your laboratory.</p>

#### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Store at room temperature.
<b>存储溶液</b>	<p>pH: 6.00</p> <p>Constituents: Citrate Buffer, 1.5% 2-butoxyethanol, 1% Hydrochloric acid</p> <p>Note: Buffer 100x concentrated, with detergent.</p>
<b>相关性</b>	Antigen Retrieval Buffer (100X Citrate Buffer pH 6.0) is a special solution that enables rehydrate and target retrieval in formalin-fixed, paraffin-embedded tissue sections in one step. It is optimal for use with primary antibodies that require citrate buffer (pH=6.0) pretreatment.

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab94674** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors