

Product datasheet

Anti-RAB23 antibody ab98164

1 图像

概述

产品名称	Anti-RAB23抗体
描述	兔多克隆抗体to RAB23
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Zebrafish 
免疫原	Synthetic peptide corresponding to a region within N-terminal amino acids 35-84 (TKDYKKTIGVDFLERQIQVNDEDVRLMLWDTAGQEEFDATKAYYRGAQ A) of Human RAB23 (NP_057361). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Fetal muscle lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab98164** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 27 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

疾病相关	Defects in RAB23 are the cause of acrocephalopolysyndactyly type 2 (ACPS2) [MIM:201000]. A syndrome characterized by craniosynostosis, polysyndactyly, obesity, and cardiac defects.
序列相似性	Belongs to the small GTPase superfamily. Rab family.
细胞定位	Cell membrane.

图片



Anti-RAB23 antibody (ab98164) at 1 µg/ml +
Fetal muscle lysate at 10 µg

Predicted band size: 27 kDa

12% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors