

Product datasheet

Anti-C14orf148 antibody ab98324

1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-C14orf148抗体 |
| 描述 | 兔多克隆抗体to C14orf148 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Human 预测可用于: Mouse, Cow, Cat, Dog |
| 免疫原 | Synthetic peptide corresponding to a region within internal sequence amino acids 179-228 (KLLLNHTNIL RPQYQYDEDS VSVWGANKGV IAALQDPTIL QATCPYSPAG) of Human C14orf148 (NP_001106946). Run BLAST with ExPASy Run BLAST with NCBI |
| 阳性对照 | HT1080 cell lysate |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: None Constituents: 2% Sucrose, PBS |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab98324** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|----|
|----|------|----|

| 应用 | Ab评论 | 说明 |
|----|------|----|
|----|------|----|


WB Use a concentration of 1 µg/ml. Predicted molecular weight: 40 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能 Probable oxidoreductase.

序列相似性 Belongs to the pyrroline-5-carboxylate reductase family.

图片



Anti-C14orf148 antibody (ab98324) at 1 µg/ml
+ HT1080 cell lysate at 10 µg

Predicted band size: 40 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors