

Product datasheet

Anti-DIP2A antibody ab98858

1 图像

概述

产品名称	Anti-DIP2A抗体
描述	兔多克隆抗体to DIP2A
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog
免疫原	Synthetic peptide corresponding to a region within the internal sequence amino acids 287-336 of Human DIP2A (NP_055966). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	HT1080 cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab98858** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 170 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能	May provide positional cues for axon pathfinding and patterning in the central nervous system.
组织特异性	Low expression in all tissues tested.
序列相似性	Belongs to the DIP2 family.
细胞定位	Nucleus.

图片



Anti-DIP2A antibody (ab98858) at 1 µg/ml (in 1 X PBS, pH 7.4 containing 5% dry milk) + HT1080 cell lysate at 10 µg

Secondary

Goat anti rabbit IgG or anti-mouse IgG -HRP at 1/10000 dilution

Developed using the ECL technique.

Predicted band size: 170 kDa

Observed band size: >170 kDa

6-18% gel

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors