

Product datasheet

Rat monoclonal [LO-MG3-7] Anti-Mouse IgG3 heavy chain (HRP) ab99646

概述

产品名称	大鼠单克隆抗体[LO-MG3-7]抗小鼠IgG3 heavy chain (HRP)
宿主	Rat
靶标种属	Mouse
特异性	ab99646 reacts with the gamma 3 heavy chain of mouse IgG3. Does not react with other immunoglobulin isotypes.
经测试应用	适用于: ICC/IF, ELISA
免疫原	Other Immunogen Type corresponding to IgG3.
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: None Constituents: 50% PBS, 50% Glycerol, pH 7.4
纯度	Purified IgM
克隆	单克隆
克隆编号	LO-MG3-7
同种型	IgM
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab99646** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		Use at an assay dependent dilution.

应用	Ab评论	说明
ELISA		1/1000 - 1/2000.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors