

Product datasheet

Mouse monoclonal [A-1] Anti-Chicken IgA heavy chain (FITC) ab99727

概述

产品名称	小鼠单克隆抗体[A-1]抗鸡IgA heavy chain (FITC)
宿主	Mouse
靶标种属	Chicken
特异性	ab99727 reacts with the alpha heavy chain of Chicken IgA.
经测试应用	适用于: ICC/IF, IHC-Fr, ELISA, Flow Cyt
偶联物	FITC. Ex: 493nm, Em: 528nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS
纯度	IgG fraction
克隆	单克隆
克隆编号	A-1
同种型	IgG2b
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab99727** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration. Acetone fixed sections
ELISA		Use at an assay dependent concentration.

应用	Ab评论	说明
Flow Cyt		Use 1µg for 10 ⁶ cells. ab18427 - Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors