

Product datasheet

Mouse monoclonal [A9604D2] Anti-Human IgA2 Fc (FITC) ab99799

概述

产品名称	小鼠单克隆抗体[A9604D2]抗人IgA2 Fc (FITC)
宿主	Mouse
靶标种属	Human
特异性	ab99799 reacts with the Fc portion of the heavy chain of IgA ₂
经测试应用	适用于: ELISA, WB, Dot blot, IHC-Fr, ICC/IF
免疫原	Other Immunogen Type corresponding to IgA2.
偶联物	FITC. Ex: 493nm, Em: 528nm

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	A9604D2
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab99799** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Predicted molecular weight: 37 kDa.

应用	Ab评论	说明
Dot blot		Use at an assay dependent dilution. for Dot and slot immunoblotting
IHC-Fr		Use at an assay dependent dilution.
ICC/IF		Use 1µg for 10 ⁶ cells.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors