

## Product datasheet

# Mouse monoclonal [HP6029] Anti-Human IgE Fc (Alkaline Phosphatase) ab99835

### 概述

产品名称	小鼠单克隆抗体[HP6029]抗人IgE Fc (Alkaline Phosphatase)
宿主	Mouse
靶标种属	Human
特异性	Reacts with the Fc portion of the heavy chain of human IgE as demonstrated by ELISA
经测试应用	适用于: ICC, Dot blot, ELISA, WB, IHC-Fr
免疫原	Other Immunogen Type corresponding to IgE.
偶联物	Alkaline Phosphatase

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: 50% Glycerol, 50mM Tris, 1mM Magnesium chloride, pH 8.0
纯度	Protein A purified
纯化说明	Source: Ascites fluid
克隆	单克隆
克隆编号	HP6029
同种型	IgG2a
轻链类型	kappa

### 应用

Our [Abpromise guarantee](#) covers the use of **ab99835** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC		Use at an assay dependent dilution.

应用	Ab评论	说明
Dot blot		Use at an assay dependent dilution.
ELISA		1/500 - 1/1000.
WB		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors